



**Australian Government**  
**Department of Health and Ageing**



Australia and New Zealand Horizon Scanning Network

**ANZHSN**

AN INITIATIVE OF THE NATIONAL, STATE AND  
TERRITORY GOVERNMENTS OF AUSTRALIA  
AND THE GOVERNMENT OF NEW ZEALAND

# **National Horizon Scanning Unit**

## **Horizon scanning prioritising summary**

**Volume 8, Number 5:**

**Vasotrac<sup>®</sup> continuous blood pressure  
monitoring.**

**February 2005**



© Commonwealth of Australia 2005

This work is copyright. You may download, display, print and reproduce this material in unaltered form only (retaining this notice) for your personal, non-commercial use or use within your organisation. Apart from any use as permitted under the Copyright Act 1968, all other rights are reserved. Requests and inquiries concerning reproduction and rights should be addressed to Commonwealth Copyright Administration, Attorney General's Department, Robert Garran Offices, National Circuit, Canberra ACT 2600 or posted at <http://www.ag.gov.au/cca>

Electronic copies can be obtained from <http://www.horizonscanning.gov.au>

Enquiries about the content of this summary should be directed to:

HealthPACT Secretariat  
Department of Health and Ageing  
MDP 106  
GPO Box 9848  
Canberra ACT 2606  
AUSTRALIA

**DISCLAIMER:** This summary is based on information available at the time of research and cannot be expected to cover any developments arising from subsequent improvements to health technologies. This summary is based on a limited literature search and is not a definitive statement on the safety, effectiveness or cost-effectiveness of the health technology covered.

The Commonwealth does not guarantee the accuracy, currency or completeness of the information in this summary. This summary is not intended to be used as medical advice and it is not intended to be used to diagnose, treat, cure or prevent any disease, nor should it be used for therapeutic purposes or as a substitute for a health professional's advice. The Commonwealth does not accept any liability for any injury, loss or damage incurred by use of or reliance on the information.

The production of this *Horizon scanning prioritising summary* was overseen by the Health Policy Advisory Committee on Technology (HealthPACT), a sub-committee of the Medical Services Advisory Committee (MSAC). HealthPACT comprises representatives from health departments in all states and territories, the Australia and New Zealand governments; MSAC and ASERNIP-S. The Australian Health Ministers' Advisory Council (AHMAC) supports HealthPACT through funding.

This *Horizon scanning prioritising summary* was prepared by Adriana Parrella from the National Horizon Scanning Unit, Adelaide Health Technology Assessment, Department of Public Health, Mail Drop 511, University of Adelaide, South Australia, 5005.

# PRIORITISING SUMMARY

**REGISTER ID:** 000144

**NAME OF TECHNOLOGY:** VASOTRAC®

**PURPOSE AND TARGET GROUP:** CONTINUOUS BLOOD PRESSURE MONITORING

**STAGE OF DEVELOPMENT (IN AUSTRALIA):**

- |   |   |
|---|---|
| <input checked="" type="checkbox"/> Yet to emerge | <input type="checkbox"/> Established  |
| <input type="checkbox"/> Experimental             | <input type="checkbox"/> Established <i>but</i> changed indication or modification of technique |
| <input type="checkbox"/> Investigational          | <input type="checkbox"/> Should be taken out of use   |
| <input type="checkbox"/> Nearly established       |   |

**AUSTRALIAN THERAPEUTIC GOODS ADMINISTRATION APPROVAL**

- |  |   |
|--|---|
| <input type="checkbox"/> Yes           | ARTG number                             |
| <input checked="" type="checkbox"/> No | <input type="checkbox"/> Not applicable |

**INTERNATIONAL UTILISATION:**

| COUNTRY       | LEVEL OF USE                 |             |                 |
|---------------|------------------------------|-------------|-----------------|
|               | Trials Underway or Completed | Limited Use | Widely Diffused |
| United States | ✓                            |             |                 |

**IMPACT SUMMARY:**

Medvave, Inc. has developed Vasotrac® to continuously monitor blood pressure. Vasotrac® was first approved for use in the United States in 2001. It is not yet available in Australia.

**BACKGROUND**




Figure 1. Vasotrac® (Printed with permission: Medvave Inc)

The Vasotrac® unit consists of a nondisposable, fluid filled cushioned sensor that is placed over the radial artery in the wrist. The sensor is connected to a monitor with a display, a control panel and processing software using a cable. The Vasotrac® can monitor blood pressure (BP) and heart rate either intermittently or continuously. When in continuous mode, the device can generate a blood pressure and heart rate measurement every 25 heartbeats and an arterial pressure waveform is displayed on the monitor and the data stored. The device

could be particularly useful for paediatric and geriatric patients and in intensive care, where obtaining BP measurements can be problematic.

Other automated, continuous measurement devices exist such as the use of an arterial catheter or an oscillometer. These techniques are more invasive and may be difficult to use on problematic patients.

#### **CLINICAL NEED AND BURDEN OF DISEASE**

It is not possible to quantify the need for blood pressure monitoring in the Australian hospital setting as measurements are routinely checked for all types of patients.

#### **DIFFUSION**

It was not possible to access information on the uptake of the Vasotrac<sup>®</sup> from the manufacturer at the time of writing this summary.

#### **COMPARATORS**

The most common method to measure blood pressure is manually utilising a sphygmomanometer. This involves the application of a cuff to the upper arm cuff that is inflated and deflated over the brachial artery to yield a numeric display of systolic and diastolic blood pressure.

An alternative for continuous monitoring is the use of an arterial catheter for beat-to-beat display of the pressure waveform and numeric display of BP. This technique is invasive and may be difficult to apply in some patient groups. There are other automated devices that take continuous measurements using an oscillometric technique with a cuff placed on an upper or lower extremity, which is automatically inflated and deflated to obtain blood pressure measurements.

#### **EFFECTIVENESS AND SAFETY ISSUES**

See complete volume of Prioritising Summaries for definitions of Levels of Evidence. The initial study (level III-3 diagnostic evidence) that assessed the effectiveness of the Vasotrac<sup>®</sup> compared it to radial artery catheterisation in four surgical centres (Belani et al 1999). Blood pressure measurements via the Vasotrac<sup>®</sup> and catheterisation of the radial artery were collected for a period of time ranging from 30 minutes to 6.2 hours, depending on the length of surgical procedure or intervention in the intensive care unit. A total of 17,468 measurements were available for analysis in 80 surgical and critically ill patients. The study reported that the Vasotrac<sup>®</sup> accurately measured BP and displayed a waveform and heart rate similar to BP measured intra-arterially. The Vasotrac<sup>®</sup> correlated highly with measurements obtained by a radial arterial catheter ( $p < 0.01$ ), systolic BP ( $r^2 = 0.93$ ) and diastolic BP ( $r^2 = 0.89$ ).

The authors reported no difficulty in measuring BP with the Vasotrac<sup>®</sup> using the radial artery, and that the device warranted further evaluation to assess performance during states of haemorrhage, hypotension, vasoconstriction and during longer periods of time than assessed in this study. The Vasotrac<sup>®</sup> was considered advantageous over the oscillometric (cuff) method as it did not interfere with pulsed oximetry monitoring equipment or the flow of intravenous fluids.

In addition, the Vasotrac<sup>®</sup> has been tested (level III-3 diagnostic evidence) on 15 patients in an Emergency Department setting (Thomas et al, 2003). Mean arterial pressures (MAP) were assessed simultaneously by the Vasotrac<sup>®</sup>, an oscillometric brachial artery cuff and an arterial line in all patients. The diagnoses for this patient group were neurological (7/15, 47%), cardiac (2/15, 13%), and other medical (6/15, 40%). There was no difference between the MAP in each of the three study groups (ANOVA  $p = 0.98$ ). The study reported that the

Vasotrac<sup>®</sup> MAP assessments demonstrated a strong statistical correlation to arterial line MAPs and performed as well as the oscillometric cuff assessment.

Finally, the Vasotrac<sup>®</sup> was tested (level III-3 diagnostic evidence) for accuracy in 15 morbidly obese patients undergoing bariatric surgery and compared to both arterial line and upper arm oscillometric cuff pressure measurement (Hager et al 2004). In this study the Vasotrac<sup>®</sup> performed as well as the arterial line and was significantly more accurate than the cuff pressure monitoring. The correlation between the arterial line and the Vasotrac<sup>®</sup> was  $r^2=0.81 \pm 0.04$  and significantly better than correlation between Vasotrac<sup>®</sup> and oscillometric cuff pressure measurement of  $r^2=0.32 \pm 0.02$  ( $p<0.001$ ).

The Vasotrac<sup>®</sup> was considered well suited to obese patients as the traditional cuff method of measuring blood pressure may be impossible to achieve or produce inaccurate measurements due to cuff size limitations.

#### **COST IMPACT**

The cost of the device in the United States is \$3,850.00 (AUD\$4,976). It was not possible to identify any information on likely cost impact of the Vasotrac<sup>®</sup> at the time of preparing this summary. Any cost impact would be the likely to occur in personnel time in placing arterial catheters, or conducting manual cuff measurements as part of regular patient observation.

#### **ETHICAL, CULTURAL OR RELIGIOUS CONSIDERATIONS**

No issues were identified/raised in the sources examined.

#### **OTHER ISSUES**

No issues were identified/raised in the sources examined.

#### **CONCLUSION:**

A non-invasive, accurate continuous blood pressure measurement device would be useful for certain patient groups and in certain medical settings. However, there is an absence of studies quantifying the impact of this technology on health or efficiency outcomes.

#### **HEALTHPACT ACTION:**

Technology will not impact significantly in terms of policy or cost burden on the Australian health system. Therefore it is recommended that this technology be archived.

#### **SOURCES OF FURTHER INFORMATION:**

Belani, K., Ozaki, M. et al (1999a). 'A new noninvasive method to measure blood pressure: results of a multicenter trial', *Anesthesiology*, 91 (3), 686-692.

Belani, K. G., Buckley, J. J. & Poliac, M. O. (1999b). 'Accuracy of radial artery blood pressure determination with the Vasotrac', *Can J Anaesth*, 46 (5 Pt 1), 488-496.

Belani, K., Komanduri, V. 2000 'Radial Artery Blood Pressure Monitoring in Children with the Vasotrac - A Preliminary Report' abstract 1278 presented American Society of Anesthesiologists Annual Meeting

Friedman, B. H., Christie, I. C. et al (2004). 'Self-reported sensitivity to continuous noninvasive blood pressure monitoring via the radial artery', *J Psychosom Res*, 57 (2), 119-121.

Hager, H., Mandadi, G., Eagon, C., Pulley, D., Kurz, A. 2004 'Intraoperative Blood Pressure Measurement on the Wrist Is More Accurate Than on the Upper Arm in Morbidly Obese Patients'. Abstract 600 presented American Society of Anesthesiologists Annual Meeting

Thomas, S., Winsor, G.R., Pang, P.S., Driscoll, K.A., Parry, B.A. 2003 'Use of a radial artery compression device for non-invasive, near-continuous blood pressure monitoring in the

emergency department.' Poster presentation, 2003 American College of Emergency Physicians Research Forum, October 12/13, Boston, MA

Medvave Inc (2004). *Blood Pressure Monitoring! Non-Invasive Direct Arterial BP Monitors by Medwave*. [Internet]. Medvave Inc. Available from: <http://www.vasotracc.com/> [Accessed 31<sup>st</sup> January 2005].

**SEARCH CRITERIA TO BE USED:**

Blood Pressure

Blood Pressure Determination

Pulse

Radial Artery/physiology