



Australian Government
Department of Health and Ageing



Australia and New Zealand Horizon Scanning Network

ANZHSN

AN INITIATIVE OF THE NATIONAL, STATE AND
TERRITORY GOVERNMENTS OF AUSTRALIA
AND THE GOVERNMENT OF NEW ZEALAND

National Horizon Scanning Unit

Horizon scanning prioritising summary

Volume 4, Number 3:

Brahms PCT Assays: For the diagnosis and control treatment of systemic bacterial infection in emergency and intensive care patients.

February 2004



© Commonwealth of Australia 2005

This work is copyright. You may download, display, print and reproduce this material in unaltered form only (retaining this notice) for your personal, non-commercial use or use within your organisation. Apart from any use as permitted under the Copyright Act 1968, all other rights are reserved. Requests and inquiries concerning reproduction and rights should be addressed to Commonwealth Copyright Administration, Attorney General's Department, Robert Garran Offices, National Circuit, Canberra ACT 2600 or posted at <http://www.ag.gov.au/cca>

Electronic copies can be obtained from <http://www.horizonscanning.gov.au>

Enquiries about the content of this summary should be directed to:

HealthPACT Secretariat
Department of Health and Ageing
MDP 106
GPO Box 9848
Canberra ACT 2606
AUSTRALIA

DISCLAIMER: This summary is based on information available at the time of research and cannot be expected to cover any developments arising from subsequent improvements to health technologies. This summary is based on a limited literature search and is not a definitive statement on the safety, effectiveness or cost-effectiveness of the health technology covered.

The Commonwealth does not guarantee the accuracy, currency or completeness of the information in this summary. This summary is not intended to be used as medical advice and it is not intended to be used to diagnose, treat, cure or prevent any disease, nor should it be used for therapeutic purposes or as a substitute for a health professional's advice. The Commonwealth does not accept any liability for any injury, loss or damage incurred by use of or reliance on the information.

The production of this *Horizon scanning prioritising summary* was overseen by the Health Policy Advisory Committee on Technology (HealthPACT), a sub-committee of the Medical Services Advisory Committee (MSAC). HealthPACT comprises representatives from health departments in all states and territories, the Australia and New Zealand governments; MSAC and ASERNIP-S. The Australian Health Ministers' Advisory Council (AHMAC) supports HealthPACT through funding.

This *Horizon scanning prioritising summary* was prepared by Adriana Parrella and Tracy Merlin from the National Horizon Scanning Unit, Adelaide Health Technology Assessment, Department of Public Health, Mail Drop 511, University of Adelaide, South Australia, 5005.

PRIORITISING SUMMARY

REGISTER ID: 0000069

NAME OF TECHNOLOGY: BRAHMS PCT ASSAYS

PURPOSE AND TARGET GROUP: DIAGNOSING AND CONTROLLING TREATMENT OF SYSTEMIC BACTERIAL INFECTION IN EMERGENCY AND INTENSIVE CARE PATIENTS

STAGE OF DEVELOPMENT (IN AUSTRALIA):

- | | |
|--|---|
| <input type="checkbox"/> Experimental | <input type="checkbox"/> Established |
| <input type="checkbox"/> Investigational | <input type="checkbox"/> Established <i>but</i> changed indication or modification of technique |
| <input checked="" type="checkbox"/> Nearly established | <input type="checkbox"/> Should be taken out of use |

AUSTRALIAN THERAPEUTIC GOODS ADMINISTRATION APPROVAL

- | | |
|------------------------------|--|
| <input type="checkbox"/> Yes | ARTG number |
| <input type="checkbox"/> No | <input checked="" type="checkbox"/> Not applicable |

INTERNATIONAL UTILISATION:

COUNTRY	LEVEL OF USE		
	Trials Underway or Completed	Limited Use	Widely Diffused
Switzerland	✓		
Australia		✓	

IMPACT SUMMARY:

Brahms have introduced a new *in vitro* assay for the quantitative determination of procalcitonin (PCT) in blood serum and plasma. PCT is a 116 amino acid protein with a sequence identical to that of the prohormone calcitonin. Under normal metabolic conditions, hormonally active calcitonin is produced and secreted in the C-cells of the thyroid gland after specific intracellular proteolytic procession of the prohormone PCT. Intact PCT is found in the blood of patients with severe bacterial infections and sepsis (Brahms, 2004).

Indications for using the assay include differential diagnosis of bacterial vs viral infection, early diagnosis of systemic bacterial and fungal infections, and monitoring for sepsis and infectious disease in high-risk patients.

The PCT assay is available in 4 formats;

- 1) **BRAHMS PCT-Q** - a rapid format for quick semi-quantitative results (available in the U.S. for investigational or research use only)
- 2) **BRAHMS PCT LIA** - a luminometer based, semi-manual, full quantitative method
- 3) **BRAHMS PCT KRYPTOR** – an automated method run on the BRAHMS Kryptor instrument (there are 7 installed in Maternal Genetic Screening Units in Australia). The KRYPTOR PCT was launched in Europe in 2002.
- 4) **LIAISON BRAHMS PCT** – an automated method run on the Diasorin Liaison instrument (not in use in Australia).

The increase in PCT in response to severe *systemic* bacterial infections may provide high diagnostic specificity. Meisner (2002) reports results from previous studies indicating that the PCT assay is 100% sensitive at 1.5 ng/ml at correctly identifying septic shock, with a specificity of 72%. The PCT assay does not react to locally limited bacterial infections, viral infections, chronic inflammatory disorders or autoimmune processes, thus enabling diagnostic differentiation between these clinical symptoms. Importantly, an increase in PCT may be detectable early in the course of severe sepsis and septic shock. The analytical sensitivity of Brahms quantitative test methods ranges between 0.06 ng/ml for Kryptor PCT and 0.1 ng/ml for PCT LIA (Meisner 2002).

A cluster randomised controlled trial (level II evidence) conducted in Switzerland assessed the effect of PCT guided diagnosis on the treatment of 234 patients admitted to hospital with suspected lower respiratory tract infections. Patients assessed with the assistance of KRYPTOR PCT assays received significantly fewer antibiotics (49% reduction, $p < 0.0001$) than patients receiving standard care (Christ-Crain et al 2004).

In Australia Brahms supplies the PCT-Q and PCT LIA formats. The Brahms PCT LIA luminometer system is currently available commercially and is the most commonly used method in Australia. It is currently not required to be listed on the Australian Register of Therapeutic Goods as the product does not diagnose a communicable disease, is not a home-use device and is not made from material of human origin. The regulations governing the listing of goods on the TGA will, however, change as of next year and the assay will be listed. The PCT LIA assay is being used in public hospitals in New South Wales, Victoria and Western Australia in emergency care and intensive care units. The quantitative tests cost approx \$12.00 – 15.00 and the semi-quantitative device cost approximately \$22.00.

In 2001-2 in Australia there were 12,688 public hospital and 2,261 private hospital separations for principal diagnosis codes A20-A49. Septicaemia is included in this range of diagnoses. This is an under-estimate of those who would actually receive the test as those found to have viral or local bacterial infections are not included. The total mortality from septicaemia in 1996 was 280 deaths in males and 315 in females (3.1 and 3.4 per 100,000).

CONCLUSION:

There is limited information currently available on the Brahms quantitative tests, however there is a small proportion of people likely to benefit from this application through earlier or more appropriate treatment.

HEALTHPACT ACTION:

It is therefore recommended that this technology be monitored.

SOURCES OF FURTHER INFORMATION:

Brahms. PCT-Molecule and Kinetics. Brahms PCT [Internet] Brahms Aktiengesellschaft, Hennigsdorf, Germany. Available from: <http://www.procalcitonin.com> [Accessed 11th March 2004].

Christ-Crain, M., Jaccard-Stolz, D. et al (2004). 'Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial', *Lancet*, 363 (9409), 600-607.

Meisner, M. (2002). 'Pathobiochemistry and clinical use of procalcitonin', *Clin Chim Acta*, 323 (1-2), 17-29.

SEARCH CRITERIA TO BE USED:

Amino Acid Sequence

Calcitonin/analysis/ biosynthesis/genetics/ metabolism

Inflammation/physiopathology

Molecular Sequence Data

Protein Precursors/analysis/ biosynthesis/genetics/ metabolism

Protein Processing, Post-Translational

Sepsis/physiopathology

Sequence Homology, Amino Acid