



**Australian Government**

**Department of Health and Ageing**



Australia and New Zealand Horizon Scanning Network

**ANZHSN**

AN INITIATIVE OF THE NATIONAL, STATE AND  
TERRITORY GOVERNMENTS OF AUSTRALIA  
AND THE GOVERNMENT OF NEW ZEALAND

## **Horizon Scanning Technology**

### **Prioritising Summary**

# **LIPOchip DNA microarray for the detection of familial hypercholesterolaemia**

**February 2009**



*Adelaide  
Health Technology  
Assessment*

© Commonwealth of Australia 2009

ISBN

Publications Approval Number:

This work is copyright. You may download, display, print and reproduce this material in unaltered form only (retaining this notice) for your personal, non-commercial use or use within your organisation. Apart from any use as permitted under the *Copyright Act 1968*, all other rights are reserved. Requests and inquiries concerning reproduction and rights should be addressed to Commonwealth Copyright Administration, Attorney General's Department, Robert Garran Offices, National Circuit, Canberra ACT 2600 or posted at <http://www.ag.gov.au/cca>

Electronic copies can be obtained from <http://www.horizonscanning.gov.au>

Enquiries about the content of the report should be directed to:

HealthPACT Secretariat  
Department of Health and Ageing  
MDP 106  
GPO Box 9848  
Canberra ACT 2606  
AUSTRALIA

**DISCLAIMER:** This report is based on information available at the time of research cannot be expected to cover any developments arising from subsequent improvements health technologies. This report is based on a limited literature search and is not a definitive statement on the safety, effectiveness or cost-effectiveness of the health technology covered.

The Commonwealth does not guarantee the accuracy, currency or completeness of the information in this report. This report is not intended to be used as medical advice and intended to be used to diagnose, treat, cure or prevent any disease, nor should it be used therapeutic purposes or as a substitute for a health professional's advice. The Commonwealth does not accept any liability for any injury, loss or damage incurred by use of or reliance the information.

The production of this Horizon scanning prioritising summary was overseen by the Health Policy Advisory Committee on Technology (HealthPACT), a sub-committee of the Medical Services Advisory Committee (MSAC). HealthPACT comprises representatives from departments in all states and territories, the Australia and New Zealand governments; and ASERNIP-S. The Australian Health Ministers' Advisory Council (AHMAC) supports HealthPACT through funding.

This Horizon scanning prioritising summary was prepared by Adrian Purins and Professor Janet Hiller from the National Horizon Scanning Unit, Adelaide Health Technology Assessment, Discipline of Public Health, School of Population Health and Clinical Practice, Mail Drop DX 650 545, University of Adelaide, Adelaide, SA, 5005.

# PRIORITISING SUMMARY

**REGISTER ID:** 000419

**NAME OF TECHNOLOGY:** LIPOCHIP

**PURPOSE AND TARGET GROUP:** GENETIC DIAGNOSIS OF FAMILIAL HYPERCHOLESTEROLAEMIA

## STAGE OF DEVELOPMENT (IN AUSTRALIA):

- |   |   |
|---|---|
| <input checked="" type="checkbox"/> Yet to emerge | <input type="checkbox"/> Established  |
| <input type="checkbox"/> Experimental             | <input type="checkbox"/> Established <i>but</i> changed indication or modification of technique |
| <input type="checkbox"/> Investigational          | <input type="checkbox"/> Should be taken out of use   |
| <input type="checkbox"/> Nearly established       |   |

## AUSTRALIAN THERAPEUTIC GOODS ADMINISTRATION APPROVAL

- |   |             |
|---|-------------|
| <input type="checkbox"/> Yes            | ARTG number |
| <input checked="" type="checkbox"/> No  |             |
| <input type="checkbox"/> Not applicable |             |

## INTERNATIONAL UTILISATION:

COUNTRY	LEVEL OF USE		
	Trials Underway or Completed	Limited Use	Widely Diffused
Spain			✓

## IMPACT SUMMARY:

Progenika-Biopharma markets the LIPOchip platform, which is designed to diagnose familial hypercholesterolaemia (FH). FH can result in premature death and is often undiagnosed until a coronary event occurs. The LIPOchip platform is a tiered system with a microarray based first tier coupled with genetic sequencing. This system is aimed at providing a definitive diagnosis (genotype) which is not possible using standard diagnosis based on clinical symptoms (phenotype). FH if diagnosed early can be treated with existing therapies.

## BACKGROUND

Familial hypercholesterolaemia is a common but under-diagnosed disorder that results in an almost 100-fold increased risk of coronary artery disease (CAD) (Bates et al 2008). FH follows an autosomal dominant pattern of inheritance, with offspring having a 50 per cent chance to inherit the disorder. Homozygotes usually have a more severe phenotype than heterozygotes, but symptom severity also results from the specific mutation(s) in each subject (Emery et al 2007).

FH is a genetic disorder resulting mainly from mutations in the low-density lipoprotein receptor (LDLR) gene that encodes the LDL receptor protein, or the apolipoprotein B (ApoB) gene. Other mutations in different genes are known but are very rare. Although over 1,000 mutations have been identified, many are rare or family specific. FH is rarely diagnosed until a coronary event occurs in a proband, which may result in other family members being screened for their FH status. The clinical symptoms of FH manifest due to the impaired uptake of cholesterol, specifically LDL, resulting in high circulating cholesterol levels and a subsequent increased risk of cardiovascular disease. Standard diagnosis of FH is performed using family history and blood lipid examination. This results in a *phenotypic* diagnosis of the patient which, due to other causes of hypercholesterolaemia, may not correlate with their genetic FH status.

The LIPOchip platform is based on a tiered mutation identification algorithm. Initially the patient is assessed using a microarray based mutation detection system. The microarray is designed to detect 203 mutations in the LDLR gene and four mutations in the ApoB gene. The subject's genetic material is extracted from a blood sample and assessed with a microarray. If the patient has any of the 207 mutations the subject is deemed FH positive. If the subject is negative their sample is subjected to a second tier of large genetic rearrangement analysis<sup>1</sup>. If this tier is also negative then a third tier consisting of sequencing of the LDLR gene to assess the subject for novel mutations. If this step is also negative then the subject is deemed FH negative.

### **CLINICAL NEED AND BURDEN OF DISEASE**

FH is a cause of premature cardiac death, with most cases not diagnosed until the primary cardiac event. It is estimated that 5-10 per cent of CAD before the age of 55 results from FH. It is thought that there may be 40,000 FH cases in Australia, with only 20 per cent correctly diagnosed and less than 10 per cent are treated with existing therapy. The prevalence of FH in the general population is thought to be 1 in 500, with sub-populations of certain ancestries having a higher prevalence, up to 1 in 70 for those of Afrikaaner origins (Emery et al 2007).

### **DIFFUSION**

The LIPOchip is not in use within Australia at the time of writing.

### **COMPARATORS**

The standard diagnostic method for FH is based on clinical testing and family history. The criteria for the standard diagnosis of FH are:<sup>2</sup>

---

<sup>1</sup> This analysis will detect if the subject's genetic material has undergone any large scale rearrangements or deletions. This type of mutation would be expected to produce a negative result in the first tier as the subject's target sequence is missing or significantly damaged.

<sup>2</sup> Adapted from Emery, J., Barlow-Stewart, K. et al (2007). 'Genetics and preventive health care', *Aust Fam Physician*, 36 (10), 808-811.

- a. Known FH DNA mutation;
- b. Tendon xanthomas in patient or first/second degree relative;
- c. Family history CHD <50 years of age in second degree relative or <60 years of age in first degree relative;
- d. Family history of cholesterol >7.5 mmol/L in first or second degree relative;
- e. Cholesterol >7.5 mmol/L (adult) or >6.7 mmol/L (age <16 years); and
- f. LDL-C >4.9 mmol/L (adult) or >4.0 mmol/L (age <16 years).

Definite FH: (e or f) + a

Probable FH: (e or f) + b

Possible FH: (e or f) + (c +d)

It is recommended that when a FH proband is diagnosed that cascade screening of relatives occurs as there is a 50 per cent chance of first degree relations being FH positive(Emery et al 2007).

#### **SAFETY AND EFFECTIVENESS ISSUES**

Several studies have used the LIPOchip to diagnose FH affected individuals. The LIPOchip has several versions with the array being updated yearly to be able to detect newly discovered mutations, and as such only the most recent studies are included.

A study using the LIPOchip investigated 825 consecutive subjects attending three Spanish lipid clinics (level IV diagnostic evidence). Subjects were suspected of having FH due to: very high familial total or LDL cholesterol levels; possible family history of premature CAD; and possible tendon xanthomas<sup>3</sup>. The version of the LIPOchip used screened for 203 mutations in the LDLR gene and four in ApoB. Subjects who had negative results were subsequently assessed for large genetic rearrangement analysis. If they were negative for this also, then genetic sequencing of the promoter, 18 exons and flanking introns was performed to uncover new mutations. Overall 459 subjects were found to have FH-causing mutations. Those who were positive for mutations were found to be more likely to have a family history of tendon xanthomas (OR 7.779 (95% CI [3.639–16.712])), more likely to have tendon xanthomas (OR 3.675 (95 % CI [1.583–8.528])), and more likely to be female (OR 1.966 (95 % CI [1.059–3.649])). No problems were reported with using the LIPOchip (Civeira et al 2008b).

A second study by the same researchers investigated the overlap of FH with familial combined hyperlipidaemia (FCH)<sup>4</sup>. The population consisted of 143 FCH patients

<sup>3</sup> Subcutaneous depositions of fat associated with tendons, especially the Achilles tendon, or those of the hands and feet. Caused by severe hypercholesterolaemia.

<sup>4</sup> The diagnostic criteria for FCH are equivocal and overlap with those of FH, and include hyperlipidaemia and hypercholesterolaemia. FCH is linked with obesity and diabetes, which is not necessarily true for the general FH positive population. The authors therefore aimed to definitively



**REFERENCES:**

Bates, T. R., Burnett, J. R. et al (2008). 'Detection of familial hypercholesterolaemia: a major treatment gap in preventative cardiology', *Heart Lung Circ*, 17 (5), 411-413.

Civeira, F., Jarauta, E. et al (2008a). 'Frequency of low-density lipoprotein receptor gene mutations in patients with a clinical diagnosis of familial combined hyperlipidemia in a clinical setting', *J Am Coll Cardiol*, 52 (19), 1546-1553.

Civeira, F., Ros, E. et al (2008b). 'Comparison of genetic versus clinical diagnosis in familial hypercholesterolemia', *Am J Cardiol*, 102 (9), 1187-1193, 1193 e1181.

Emery, J., Barlow-Stewart, K. et al (2007). 'Genetics and preventive health care', *Aust Fam Physician*, 36 (10), 808-811.

**SEARCH CRITERIA TO BE USED:**

DNA Mutational Analysis

Receptors, LDL/ blood

Cholesterol, LDL/ blood

Receptors, Lipoprotein/ blood