



Australian Government

Department of Health and Ageing



Australia and New Zealand Horizon Scanning Network

ANZHSN

AN INITIATIVE OF THE NATIONAL, STATE AND
TERRITORY GOVERNMENTS OF AUSTRALIA
AND THE GOVERNMENT OF NEW ZEALAND

Horizon Scanning Technology Prioritising Summary

EasyOne™ spirometer for the diagnosis and management of chronic respiratory disease and asthma

May 2007



© Commonwealth of Australia 2007

[add ISSN]

[add Publications Approval Number]

This work is copyright. You may download, display, print and reproduce this material in unaltered form only (retaining this notice) for your personal, non-commercial use or use within your organisation. Apart from any use as permitted under the Copyright Act 1968, all other rights are reserved. Requests and inquiries concerning reproduction and rights should be addressed to Commonwealth Copyright Administration, Attorney General's Department, Robert Garran Offices, National Circuit, Canberra ACT 2600 or posted at <http://www.ag.gov.au/cca>

Electronic copies can be obtained from <http://www.horizonscanning.gov.au>

Enquiries about the content of this summary should be directed to:

HealthPACT Secretariat
Department of Health and Ageing
MDP 106
GPO Box 9848
Canberra ACT 2606
AUSTRALIA

DISCLAIMER: This summary is based on information available at the time of research and cannot be expected to cover any developments arising from subsequent improvements to health technologies. This summary is based on a limited literature search and is not a definitive statement on the safety, effectiveness or cost-effectiveness of the health technology covered.

The Commonwealth does not guarantee the accuracy, currency or completeness of the information in this summary. This summary is not intended to be used as medical advice and it is not intended to be used to diagnose, treat, cure or prevent any disease, nor should it be used for therapeutic purposes or as a substitute for a health professional's advice. The Commonwealth does not accept any liability for any injury, loss or damage incurred by use of or reliance on the information.

The production of this *Horizon scanning prioritising summary* was overseen by the Health Policy Advisory Committee on Technology (HealthPACT), a sub-committee of the Medical Services Advisory Committee (MSAC). HealthPACT comprises representatives from health departments in all states and territories, the Australia and New Zealand governments; MSAC and ASERNIP-S. The Australian Health Ministers' Advisory Council (AHMAC) supports HealthPACT through funding.

This *Horizon scanning prioritising summary* was prepared by Linda Mundy and Janet Hiller from the National Horizon Scanning Unit, Adelaide Health Technology Assessment, Discipline of Public Health, Mail Drop 511, University of Adelaide, South Australia, 5005.

PRIORITISING SUMMARY

REGISTER ID: 000302

NAME OF TECHNOLOGY: EASYONE™ SPIROMETER

PURPOSE AND TARGET GROUP: FOR THE DIAGNOSIS AND MANAGEMENT OF CHRONIC RESPIRATORY DISEASE AND ASTHMA

STAGE OF DEVELOPMENT (IN AUSTRALIA):

- | | |
|--|---|
| <input type="checkbox"/> Yet to emerge | <input type="checkbox"/> Established |
| <input type="checkbox"/> Experimental | <input type="checkbox"/> Established <i>but</i> changed indication or modification of technique |
| <input type="checkbox"/> Investigational | <input type="checkbox"/> Should be taken out of use |
| <input checked="" type="checkbox"/> Nearly established | |

AUSTRALIAN THERAPEUTIC GOODS ADMINISTRATION APPROVAL

- | | | |
|---|-------------|-------|
| <input checked="" type="checkbox"/> Yes | ARTG number | 81644 |
| <input type="checkbox"/> No | | |
| <input type="checkbox"/> Not applicable | | |

INTERNATIONAL UTILISATION:

COUNTRY	LEVEL OF USE		
	Trials Underway or Completed	Limited Use	Widely Diffused
Australia		✓	
New Zealand		✓	
Multi-centre – South America	✓		
United States	✓		

IMPACT SUMMARY:

The EasyOne™ spirometer is manufactured by ndd Medizintechnik (Switzerland) and distributed in Australia by Niche Medical. The EasyOne™ spirometer is designed for the diagnosis and management of chronic respiratory disease and asthma. The EasyOne™ is registered on the TGA and received approval from the FDA in 2000.

BACKGROUND

Undiagnosed airway obstruction is common in the general population. In the adult population not diagnosed with asthma, the main cause (90%) of airflow limitation is chronic obstructive pulmonary disease (COPD) due to cigarette smoking. COPD may also be caused by airborne irritants such as coal dust, air pollution from burning wood, asbestos or solvents (Enright & Kaminsky 2003). The main symptoms of COPD include

dyspnea (shortness of breath) which may be accompanied by wheezing, persistent cough with sputum and chronic bronchitis. A rare complication may be the partial or complete collapse of the lung, usually due to an obstruction of a bronchus (Wikipedia 2007). COPD is a progressive condition if left untreated or exposure to the causative agent is not halted, however COPD is preventable and treatable but not fully reversible (GOLD executive committee 2006).

The gold standard for the detection of airflow limitation and the diagnosis of COPD in its preclinical stage is spirometry. Spirometry is a non-invasive test which measures the volume of exhaled air, after a maximal inhalation, during a specified period of time (Enright & Kaminsky 2003). The results of a spirometry test are expressed as the FEV₁/FVC ratio, where FEV₁ is the forced expiratory volume in one second and FVC is the forced vital capacity, or the total amount of air forcibly exhaled after full inspiration. The severity of COPD is classified according to the following ratios:

Stage I (mild):	FEV ₁ /FVC <0.70, FEV ₁ ≥80% predicted ¹ ;
Stage II (moderate):	FEV ₁ /FVC <0.70, 50% ≤ FEV ₁ <80% predicted;
Stage III (severe):	FEV ₁ /FVC <0.70, 30% ≤ FEV ₁ <50% predicted; and
Stage IV (very severe)	FEV ₁ /FVC <0.70, FEV ₁ <30% predicted or FEV ₁ <50% predicted plus chronic respiratory failure (GOLD executive committee 2006).

Spirometry may be conducted by general practitioners in their practices using hand held spirometers. Recent updated international guidelines have recommended that the accuracy of diagnostic spirometers be checked (ie calibrated) at least once daily using a syringes with a volume of at least 3 litres. A recent survey found that although 64 per cent of general practices in Australia owned a spirometer, only 22 per cent of these practices use a calibration syringe to calibrate the spirometer and that only 1.5 per cent actually perform this task daily. Failure to calibrate may produce spurious results (Walters et al 2006).



The EasyOne™ is a handheld spirometer which utilises an ultrasonic sensor to measure airflow (Figure 1). The EasyOne™ has no moving parts, therefore its accuracy is not dependent on mechanical function and it does not require calibration. In addition, the EasyOne™ is not dependent on the measurement of variables such as pressure or displacement volume (Niche Medical 2003; Walters et al 2006).

Figure 1 The EasyOne™
(printed with permission Niche Medical)

¹ Predicted values are based on values from normal population and are calculated using the patient's height, weight and gender

Disposable spirettes (mouthpiece) with a bacterial shield are fitted for each new patient. Transducers are fitted on both sides of the spirette cavity which receive sound in alternating directions (Figure 2). When gas flow (breath sample) is present in the tube, a pulse that is travelling with the flow (downstream) (Figure 2b) reaches the opposite transducer faster than the pulse that is travelling against the flow (upstream), which is slower and takes a longer time to reach the opposite transducer (Figure 2c). The transit time of the sound pulses are measured and the gas flow through the spirette is calculated from these times. The EasyOne™ can be used as a stand alone device which is capable of storing up to 700 patient results, or the results may be uploaded via a USB port for data storage (Niche Medical 2003).




Figure 2 Gas flow into the EasyOne™ (Niche Medical 2003)

CLINICAL NEED AND BURDEN OF DISEASE

Chronic respiratory diseases are a major health problem in Australia, in 2001 it was estimated 5.8 million Australians suffered at least one long-term respiratory condition. The two major manifestations of chronic respiratory disease are asthma and COPD.

COPD is a group of long-term lung diseases, the major forms of which are emphysema and chronic bronchitis, characterized by progressive shortness of breath associated with abnormal inflammation of the lungs. The major cause of COPD is smoking; other causes are exposure to environmental hazards and genetic predisposition. In Australia COPD is a major cause of mortality and morbidity, being the underlying cause of death in 5,378 deaths in 2003 and resulting in more than 53,000 hospital separations averaging 7.5 days per separation (AIHW 2005b). Although the prevalence of COPD is difficult to determine, in 2001 it was estimated that 3.5 percent of Australians suffered from emphysema or bronchitis, although this is likely to be an underestimate as, generally, only

advanced COPD is diagnosed. Females are affected slightly more than males and increasing age is also a risk factor for diagnosis of COPD.

Asthma is a chronic disease characterised by inflammation of the airways, resulting in constriction of the airways leading to chest tightness, wheezing and shortness of breath. It is believed this cascade of events is the result of exposure to specific triggers, such as allergens, exertion, stress etc. Since the mid 1990s the prevalence of asthma in Australia may have plateaued, despite this, Australia's rate remains one of the highest levels in the world at 10-12% in children and 14-16% in adults (AIHW 2005a). In 2002-03 asthma was the principal diagnosis in 37,230 hospital separations averaging 2.5 days per separation (AIHW 2005b). Asthma was the underlying cause of death in 314 (1.6 per 100,000 population, 95% CI [1.4, 1.8]) and an associated cause of death in 934 cases in 2003 (AIHW 2005a).

Together, asthma and COPD are estimated to cost the Australian health care system over one billion dollars a year (AIHW 2005b).

DIFFUSION

The EasyOne™ is in limited use in several hospitals and clinics in Australia and New Zealand. Applications of the EasyOne™ spirometer in Australia now include: COPD, asthma and epidemiological research, use in emergency departments, respiratory wards, respiratory investigation units, COPD outreach programmes, respiratory outpatient clinics, cystic fibrosis clinics, occupational health clinics, community health centres and general practice.

COMPARATORS

Spirometers other than the EasyOne™ may be used to diagnose COPD, however these require daily calibration. Other methods for measuring airflow limitation include a simple forced expiratory time using a stethoscope and timing the expiration. If the forced expiratory time exceeds six seconds it is likely that the patient has airflow limitation, however this method results in a high misclassification rate. Peak expiratory flow (PEF) is low in patients with airflow limitation. Clinical practice guidelines recommend PEF to be used to *manage* but not to *diagnose* asthma. PEF is insensitive to airflow limitation, and intra-subject and inter-subject variability is approximately twice that of FEV. Airway obstruction increases airway resistance, which can be measured using plethysmography, a more expensive technique than spirometry and may also give variable results (Enright & Kaminsky 2003).

EFFECTIVENESS AND SAFETY ISSUES

Walters et al (2006) examined the reproducibility and accuracy of results obtained with six EasyOne™ spirometers used in general practice for between 15.1 and 26.6 weeks (mean time 23.9 weeks) (level IV diagnostic evidence). Calibration checks were carried out on the spirometers with a 3-litre calibration syringe with a certified accuracy of ± 0.5

per cent. A total of 1,041 patients were tested using the six spirometers and 75 paired syringes calibrations were performed using both a single dedicated calibration spirette and new spirette randomly chosen. All calibrations met the accuracy criteria for both expiration and inspiration checks (Table 1), however, the mean volumes for both expiration and inspiration were significantly higher ($p < 0.001$) when using a random spirette when compared to the dedicated spirette. The authors conclude that this study supports the manufacturers claim that the EasyOne™ does not require calibration (Walters et al 2006).

Table 1 Results of paired calibration checks with dedicated and random spirettes

Dedicated spirette		Random spirette	
Volume deviation from 3-L		Volume deviation from 3-L	
Absolute (mean ± SD)	% (mean ± SD)	Absolute (mean ± SD)	% (mean ± SD)
Expiration	Expiration	Expiration	Expiration
0.011 ± 0.033	0.373 ± 1.11	0.046 ± 0.034	1.523 ± 1.125
Inspiration	Inspiration	Inspiration	Inspiration
-0.036 ± 0.03	-1.209 ± 0.992	0.003 ± 0.039	0.09 ± 1.303

SD = standard deviation

A similar study conducted in five Latin American cities collected calibration data from 70 EasyOne™ spirometers used for a population based survey of COPD (level IV diagnostic evidence). This study reported that 97 per cent of the calibration volumes were within ± 64ml (2.1%) of the 3-litre calibration signal, and that 98% were within ± 50 ml (1.7%) (Perez-Padilla et al 2006).

A study by Liistro et al (2006) compared the ease of use and accuracy of 10 spirometers, including the EasyOne™, which could be routinely used in general practice. These spirometers were compared to standard diagnostic spirometers the Vmax 20C, a flow sensing spirometer, and the Morgan TLC, a volume sensing device (level III-2 diagnostic evidence). Devices were calibrated daily with a 3-litre syringe. All testing was repeated in three laboratories following the same protocol. Pulmonary function tests (PFTs) were performed in healthy non-smoking subjects and COPD patients, with each performing five successive forced expirations, resulting in a within subject standard deviation (Sw). The larger the Sw, the lower the precision of the spirometer. A total of 399 different subjects were studied in the three centres (Liistro et al 2006).

The 95% limits of precision ($1.96 \times Sw$) for forced vital capacity (FVC) in healthy subjects for the EasyOne™ was >200ml at 0.25, but was 0.19 for the forced expiratory volume (FEV), compared to the standard spirometer Sw values of 0.15 and 0.14 respectively. Only two other spirometers had a FVC Sw value higher than that reported for the EasyOne™ (0.34 and 0.43), with the remaining seven spirometers reporting a mean of 0.17 ± 0.016 . One spirometer reported a FEV Sw value of 0.29 (limits of precision >200ml) with the remaining nine spirometers (including the EasyOne™) reporting a mean of 0.15 ± 0.03 . The EasyOne™ did not return a significant bias from

zero (when compared to standard spirometers) for FVC (0.07 95%CI [0.00, 0.15], limits of agreement² -0.43 to 0.57, $r = -0.15$), or FEV₁/FVC (1.8 95%CI [0.6, 3.0], limits of agreement -6.4 to 10.0, $r = -0.26$). The EasyOne™ did present with a significant bias ($p < 0.01$) for FEV₁, however this remained in the acceptable limit of ≤ 100 ml. Note: only one other device presented a significant bias for FVC, four devices for FEV₁ and two devices for FEV₁/FVC. In conclusion, the repeatability of FVC and FEV₁/FVC measurements was good for the EasyOne™, was poorer for FEV₁ measurements but the level of agreement between the EasyOne™ and standard devices was good. The EasyOne™ was found not to be as precise measuring FVC but precision was good measuring FEV₁, when compared to standard spirometers (Liistro et al 2006).

A study by Schoh et al (2002) conducted at a health fair in a community hospital over two days, compared the spirometry results obtained with the EasyOne™ to a standard spirometer, the Vmax. A minimum of three tests were conducted with both the EasyOne™ and the Vmax devices. A total of 394 subjects were tested over the two days with the EasyOne™ and although confirmatory testing with the Vmax was offered to all participants, only those with an airflow limitation tended to take up this option. As a result, only 115 participants had a spirometry assessment with both the EasyOne™ and the Vmax. Correlations between the EasyOne™ and the Vmax were good for FEV₁ ($r = 0.93$), and for FEV₆³ ($r = 0.96$) and FEV₁/FEV₆ ($r = 0.72$) ($p = 0.001$ for all comparisons). The 95% limits of agreement were -0.18 to 0.69 for FEV₁, -0.24 to 0.81 for FEV₆, and -0.12 to 0.13 for FEV₁/FEV₆. The EasyOne™ device was found to slightly underestimate the results of the Vmax which the authors felt was a result of a learning curve effect as the spirometry technicians were well versed in the use of the standard Vmax device (level III-2 diagnostic evidence) (Schoh et al 2002).

COST IMPACT

The EasyOne™ spirometer costs approximately \$3,000 (Walters et al 2006).

ETHICAL, CULTURAL OR RELIGIOUS CONSIDERATIONS

No issues were identified/raised in the sources examined.

OTHER ISSUES

No issues were identified/raised in the sources examined.

CONCLUSION:

Based on the available evidence it appears that the EasyOne™ is a reliable device which is easy to use in a general practitioner setting. Readings from the EasyOne™ correlate well with readings from standard spirometers and the EasyOne™ device has the added

² 95% limits of agreement = mean difference between the two spirometers ± 1.96 standard deviations

³ FEV₆ = forced expiratory volume in the first 6 seconds

