



Australian Government

Department of Health and Ageing



Australia and New Zealand Horizon Scanning Network

ANZHSN

AN INITIATIVE OF THE NATIONAL, STATE AND
TERRITORY GOVERNMENTS OF AUSTRALIA
AND THE GOVERNMENT OF NEW ZEALAND

National Horizon Scanning Unit

Horizon scanning prioritising summary

Volume 5, Number 1:

**I-STAT[®] Cardiac Troponin I (cTnl) test:
for the assessment of biomarkers for acute
myocardial infarction in patients presenting
to emergency departments.**

May 2004



© Commonwealth of Australia 2005

[add ISSN]

[add Publications Approval Number]

This work is copyright. You may download, display, print and reproduce this material in unaltered form only (retaining this notice) for your personal, non-commercial use or use within your organisation. Apart from any use as permitted under the Copyright Act 1968, all other rights are reserved. Requests and inquiries concerning reproduction and rights should be addressed to Commonwealth Copyright Administration, Attorney General's Department, Robert Garran Offices, National Circuit, Canberra ACT 2600 or posted at <http://www.ag.gov.au/cca>

Electronic copies can be obtained from <http://www.horizonscanning.gov.au>

Enquiries about the content of this summary should be directed to:

HealthPACT Secretariat
Department of Health and Ageing
MDP 106
GPO Box 9848
Canberra ACT 2606
AUSTRALIA

DISCLAIMER: This summary is based on information available at the time of research and cannot be expected to cover any developments arising from subsequent improvements to health technologies. This summary is based on a limited literature search and is not a definitive statement on the safety, effectiveness or cost-effectiveness of the health technology covered.

The Commonwealth does not guarantee the accuracy, currency or completeness of the information in this summary. This summary is not intended to be used as medical advice and it is not intended to be used to diagnose, treat, cure or prevent any disease, nor should it be used for therapeutic purposes or as a substitute for a health professional's advice. The Commonwealth does not accept any liability for any injury, loss or damage incurred by use of or reliance on the information.

The production of this *Horizon scanning prioritising summary* was overseen by the Health Policy Advisory Committee on Technology (HealthPACT), a sub-committee of the Medical Services Advisory Committee (MSAC). HealthPACT comprises representatives from health departments in all states and territories, the Australia and New Zealand governments; MSAC and ASERNIP-S. The Australian Health Ministers' Advisory Council (AHMAC) supports HealthPACT through funding.

This *Horizon scanning prioritising summary* was prepared by Adriana Parrella and Tracy Merlin from the National Horizon Scanning Unit, Adelaide Health Technology Assessment, Department of Public Health, Mail Drop 511, University of Adelaide, South Australia, 5005.

PRIORITISING SUMMARY

REGISTER ID: 000092

NAME OF TECHNOLOGY: I-STAT[®] CARDIAC TROPONIN I (cTNI) TEST

PURPOSE AND TARGET GROUP: ASSESSMENT OF BIOMARKERS FOR ACUTE MYOCARDIAL INFARCTION FOR PATIENTS PRESENTING TO EMERGENCY DEPARTMENTS

STAGE OF DEVELOPMENT (IN AUSTRALIA AND/OR NEW ZEALAND):

- | | |
|--|---|
| <input type="checkbox"/> Yet to emerge | <input type="checkbox"/> Established |
| <input type="checkbox"/> Experimental | <input type="checkbox"/> Established <i>but</i> changed indication or modification of technique |
| <input type="checkbox"/> Investigational | <input type="checkbox"/> Should be taken out of use |
| <input checked="" type="checkbox"/> Nearly established | |

AUSTRALIAN THERAPEUTIC GOODS ADMINISTRATION APPROVAL

- | | | |
|---|---|-------|
| <input checked="" type="checkbox"/> Yes | ARTG number | 73566 |
| <input type="checkbox"/> No | <input type="checkbox"/> Not applicable | |

The I-STAT System Analyzer is available in over 400 sites in Australia for other blood testing applications.

INTERNATIONAL UTILISATION:

COUNTRY	LEVEL OF USE		
	Trials Underway or Completed	Limited Use	Widely Diffused
United States	✓		

IMPACT SUMMARY:

The I-STAT Corporation provides the I-STAT[®] Cardiac Troponin I (cTNI) test with the aim of assessing Troponin I levels (in heparinised whole blood or plasma samples) in patients with chest pain presenting to Emergency Departments. The manufacturer states it can be used in the diagnosis and monitoring of acute myocardial infarction and in assessing the risk of mortality.

The I-STAT[®] test is contained in a single test cartridge and is used with the I-STAT System Analyzer. The cartridge containing the blood sample is inserted into the I-STAT Analyzer and levels are measured automatically. The I-STAT Analyzer is not a new device: it is currently used for other blood testing.

BACKGROUND

Elevations of blood serum Troponin T and I values may be used to rule out acute myocardial infarction (AMI) before discharging patients from emergency departments.

A systematic review of Troponin T and I levels for diagnosing AMI in Emergency Departments included eleven level I criteria (prospective studies of patients, diagnosed with WHO criteria for AMI, and assessment blinded to Troponin results) and an additional 8 level II criteria (as level I but nonconsecutive or unspecified) studies. Diagnostic sensitivity and specificity for different points in time from patient arrival in Emergency Department or from onset of pain are presented in this review (Ebell et al. 2000). The review concluded that sensitivity to Troponin T and I increases from 10% to 45% within one hour of the onset of pain to more than 90% at eight or more hours. Specificity declines gradually from 87% to 80% from 1 to 12 hours after the onset of chest pain for Troponin T and is approximately 95% for Troponin I after 1 to 12 hours. The authors conclude that the use of Troponin T and I as biomarkers of AMI needs to be interpreted according to the number of hours from onset of pain. Currently, manual testing of Troponin I biomarkers takes over one hour, the advantage of the I-STAT[®] Cardiac Troponin test is that it is performed at point of care and produces a result in 10 minutes.

CLINICAL NEED AND BURDEN OF DISEASE

In Australia 7,484 males and 6,959 females died from AMI (777.1 and 711 per million people respectively) in the year 2001 (AIHW 2004). There were 40,338 hospital separations by principal diagnosis in ICD-10-AM for AMI in 2001-02.

DIFFUSION

Troponin testing for AMI is currently conducted in hospital laboratories in Australia (Davey, 2003). The I-STAT (cTnI) test cartridges were recently introduced to Australia in March 2004 (personal communication, I-STAT corp. company spokesperson). Given that Troponin testing is standard in Australian hospitals for diagnosing AMI, the increased speed of acquiring results through an I-STAT Analyzer compared to standard laboratory testing, suggests that the uptake of the I-STAT (cTnI) test cartridges will be quick.

COMPARATORS

Until recently the World Health Organisation criteria for diagnosing AMI included elevation of creatine kinase (CK), myocardial bound creatine kinase (CK – MB) blood levels, along with electrocardiographic changes and a clinical history compatible with ischaemia. The measurement of total CK is no longer recommended for the routine *diagnosis* of AMI, because of the wide tissue distribution of this enzyme.

The protocols for AMI were revised and consensus guidelines from the American Heart Association, American College of Cardiology and European Society of Cardiology have endorsed the use of serum or plasma cardiac Troponins for the diagnosis of AMI (Antman et al 2000, Jaffe et al. 2000).

The revision of guidelines for biomarkers in diagnosing AMI states that cardiac Troponin (I or T) is the preferred biomarker for myocardial damage due to its nearly absolute specificity, as well as high sensitivity, thereby reflecting even microscopic zones of necrosis of myocardial tissue. The protocols suggest that negative Troponin levels measured at 6 hours after symptom onset and again at 12 hours can safely determine whether a patient's symptoms are caused by an AMI (Antman 2000, Wu 1999).

The I-STAT test is the first point of care test that provides rapid measurement of Troponin I blood levels in 10 minutes. The standard method of measuring Troponin I in patients presenting in Emergency departments is to analyse blood samples in a hospital laboratory: this method can take up to and over one hour.

COST IMPACT

The current costs of pathology testing for creatine kinase, Troponin or myoglobin testing in plasma or serum for MBS item numbers 66518 and 66519, are \$20.40 and \$40.85. The cost of the I-STAT cartridge is approximately \$15.00.

EFFECTIVENESS AND SAFETY

There were no studies available that assessed the effectiveness or safety of the I-STAT[®] test. Information from the American Food and Drug Administration approval process refers to clinical testing undertaken at three clinical sites comparing the I-STAT cTnl test with another Troponin test, the Stratus CS. This reports that the lower limit of detection for the I-STAT method is slightly lower at 0.02ng/ml versus 0.03ng.ml for the Stratus method. There appear to be no reports of the I-STAT cTNL compared to other Troponin tests.

ETHICAL, CULTURAL OR RELIGIOUS CONSIDERATIONS

No issues were identified/raised in the sources examined.

CONCLUSION:

There is a lack of evidence on the use of the I-STAT cTNL. In addition, it is considered that this technology doesn't differ sufficiently from existing technologies already in use.

HEALTHPACT ACTION:

Therefore it is recommended that this technology be archived.

SOURCES OF FURTHER INFORMATION:

Adams, J. E., 3rd, Schechtman, K. B. et al (1994). 'Comparable detection of acute myocardial infarction by creatine kinase MB isoenzyme and cardiac Troponin I', *Clin Chem*, 40 (7 Pt 1), 1291-1295.

Antman, E., Bassand, J.-P. et al (2000). 'Myocardial infarction redefined--a consensus document of The Joint European Society of Cardiology/American College of Cardiology committee for the redefinition of myocardial infarction: The Joint European Society of Cardiology/ American College of Cardiology Committee*', *Journal of the American College of Cardiology*, 36 (3), 959-969.

Antman, E. M., Tanasijevic, M. J. et al (1996). 'Cardiac-specific Troponin I levels to predict the risk of mortality in patients with acute coronary syndromes', *N Engl J Med*, 335 (18), 1342-1349.

Apple, F. S., Quist, H. E. & Murakami, M. M. (2004). 'Diagnostic and prognostic value of cardiac Troponin I assays in patients admitted with symptoms suggestive of acute coronary syndrome', *Arch Pathol Lab Med*, 128 (4), 430-434.

Ebell, M. H., Flewelling, D. & Flynn, C. A. (2000). 'A systematic review of Troponin T and I for diagnosing acute myocardial infarction', *J Fam Pract*, 49 (6), 550-556.

Hamm, C. W., Goldmann, B. U. et al (1997). 'Emergency room triage of patients with acute chest pain by means of rapid testing for cardiac Troponin T or Troponin I', *N Engl J Med*, 337 (23), 1648-1653.

Jaffe, A. S., Ravkilde, J. et al (2000). 'It's time for a change to a Troponin standard', *Circulation*, 102 (11), 1216-1220.

Lin, J. C., Apple, F. S. et al (2004). 'Rates of positive cardiac Troponin I and creatine kinase MB mass among patients hospitalized for suspected acute coronary syndromes', *Clin Chem*, 50 (2), 333-338.

Meyer, T., Binder, L. et al (1998). 'Superiority of combined CK-MB and Troponin I measurements for the early risk stratification of unselected patients presenting with acute chest pain', *Cardiology*, 90 (4), 286-294.

Penttila, I., Penttila, K. & Rantanen, T. (2000). 'Laboratory diagnosis of patients with acute chest pain', *Clin Chem Lab Med*, 38 (3), 187-197.

Polanczyk, C. A., Lee, T. H. et al (1998). 'Cardiac Troponin I as a predictor of major cardiac events in emergency department patients with acute chest pain', *J Am Coll Cardiol*, 32 (1), 8-14.

Wang, L., Wu, Y. & Li, L. (2000). '[Evaluation on usage of serum Troponin-I, CK-MB mass and myoglobin measurements in diagnosing acute myocardial infarction]', *Hua Xi Yi Ke Da Xue Xue Bao*, 31 (2), 239-241, 245.

Wu, A. H., Apple, F. S. et al (1999). 'National Academy of Clinical Biochemistry Standards of Laboratory Practice: recommendations for the use of cardiac markers in coronary artery diseases', *Clin Chem*, 45 (7), 1104-1121.

Zarich, S. W., Bradley, K. et al (2004). 'Minor elevations in Troponin T values enhance risk assessment in emergency department patients with suspected myocardial ischemia: analysis of novel Troponin T cut-off values', *Clin Chim Acta*, 343 (1-2), 223-229.

SEARCH CRITERIA TO BE USED:

Biological Markers/blood

Chest Pain/diagnosis

Coronary Disease/blood/ diagnosis

Myocardial Infarction/ diagnosis

Troponin I/ blood