



Australian Government
Department of Health and Ageing



Australia and New Zealand Horizon Scanning Network

ANZHSN

AN INITIATIVE OF THE NATIONAL, STATE AND
TERRITORY GOVERNMENTS OF AUSTRALIA
AND THE GOVERNMENT OF NEW ZEALAND

National Horizon Scanning Unit

Horizon scanning prioritising summary

Volume 4, Number 5:

**Uni-GoldTM Recombigen[®] HIV test: To
detect antibodies to HIV-1 in the plasma
and whole blood of individuals at risk of
HIV infection.**

February 2004



© Commonwealth of Australia 2005

This work is copyright. You may download, display, print and reproduce this material in unaltered form only (retaining this notice) for your personal, non-commercial use or use within your organisation. Apart from any use as permitted under the Copyright Act 1968, all other rights are reserved. Requests and inquiries concerning reproduction and rights should be addressed to Commonwealth Copyright Administration, Attorney General's Department, Robert Garran Offices, National Circuit, Canberra ACT 2600 or posted at <http://www.ag.gov.au/cca>

Electronic copies can be obtained from <http://www.horizonscanning.gov.au>

Enquiries about the content of this summary should be directed to:

HealthPACT Secretariat
Department of Health and Ageing
MDP 106
GPO Box 9848
Canberra ACT 2606
AUSTRALIA

DISCLAIMER: This summary is based on information available at the time of research and cannot be expected to cover any developments arising from subsequent improvements to health technologies. This summary is based on a limited literature search and is not a definitive statement on the safety, effectiveness or cost-effectiveness of the health technology covered.

The Commonwealth does not guarantee the accuracy, currency or completeness of the information in this summary. This summary is not intended to be used as medical advice and it is not intended to be used to diagnose, treat, cure or prevent any disease, nor should it be used for therapeutic purposes or as a substitute for a health professional's advice. The Commonwealth does not accept any liability for any injury, loss or damage incurred by use of or reliance on the information.

The production of this *Horizon scanning prioritising summary* was overseen by the Health Policy Advisory Committee on Technology (HealthPACT), a sub-committee of the Medical Services Advisory Committee (MSAC). HealthPACT comprises representatives from health departments in all states and territories, the Australia and New Zealand governments; MSAC and ASERNIP-S. The Australian Health Ministers' Advisory Council (AHMAC) supports HealthPACT through funding.

This *Horizon scanning prioritising summary* was prepared by Adriana Parrella and Tracy Merlin from the National Horizon Scanning Unit, Adelaide Health Technology Assessment, Department of Public Health, Mail Drop 511, University of Adelaide, South Australia, 5005.

PRIORITISING SUMMARY

REGISTER ID: 0000078

NAME OF TECHNOLOGY: UNI-GOLD™ RECOMBIGEN® HIV TEST

PURPOSE AND TARGET GROUP: TO DETECT ANTIBODIES TO HIV-1 IN THE PLASMA AND WHOLE BLOOD OF INDIVIDUALS AT RISK OF HIV INFECTION

STAGE OF DEVELOPMENT (IN AUSTRALIA):

- | | |
|---|---|
| <input checked="" type="checkbox"/> Yet to emerge | <input type="checkbox"/> Established |
| <input type="checkbox"/> Experimental | <input type="checkbox"/> Established <i>but</i> changed indication or modification of technique |
| <input type="checkbox"/> Investigational | <input type="checkbox"/> Should be taken out of use |
| <input type="checkbox"/> Nearly established | |

AUSTRALIAN THERAPEUTIC GOODS ADMINISTRATION APPROVAL

- | | |
|--|---|
| <input type="checkbox"/> Yes | ARTG number |
| <input checked="" type="checkbox"/> No | <input type="checkbox"/> Not applicable |

INTERNATIONAL UTILISATION:

COUNTRY	LEVEL OF USE		
	Trials Underway or Completed	Limited Use	Widely Diffused
USA	✓		

IMPACT SUMMARY:

Trinity Biotech provides the Uni-Gold™ Recombigen® HIV test for the detection of antibodies to HIV-1 in plasma, serum and whole blood. The test was approved in the U.S.A. in December 2003. The technology has not yet emerged in Australia.

Alternate procedures to test for HIV-1 infection are to run assays for the various virus components or by tests that assess whether an individual's immune system has produced an HIV-1-specific response. Samples are screened with an ELISA and positives are confirmed with Western Blot or Immuno-fluorescence assay. The standard laboratory test algorithm may take 48 hours to one week whereas the Uni-Gold Recombigen HIV provides results in 10 minutes (Trinity Biotech 2003). The advantage of rapid HIV testing may be demonstrated in urgent medical circumstances and in settings where individuals do not return for HIV test results.

Uni-Gold™ Recombigen® HIV detects antibodies to HIV-1 using colloidal gold to give a result. Clinical studies (level 3b evidence, see Appendix A) of 1,032 HIV-1 positive samples undertaken by Trinity Biotech ascertained sensitivity at 100% [95% C.I. 99.5%, 100%]. 1,000 of the samples were collected from individuals known to be HIV-1 sero-positive, and previously confirmed by western blot. A further 32 samples were collected from individuals from high risk populations of unknown HIV serostatus who were subsequently found to be repeatedly reactive using ELISA and Western Blot. The Uni-Gold™ Recombigen® HIV test was reactive for all samples using serum, plasma and whole blood portions (Trinity Biotech 2003).

The specificity of the Uni-Gold™ Recombigen® HIV was evaluated using fresh serum, plasma and whole blood samples. A total of 1,968 HIV-1 ELISA negative individual samples were run as serum, plasma and whole blood on Uni-Gold™ Recombigen® HIV (level 3b diagnostic evidence). 1,000 were collected from individuals of unknown HIV status in a low risk population and subsequently confirmed as negative by ELISA. Of these samples, two were reactive in initial test by plasma and serum and three by whole blood when tested by Uni-Gold™ Recombigen® HIV. In a low risk population the specificity of Uni-Gold™ Recombigen® HIV in these studies was 99.8% [95% C.I. 99.3%, 100%] for serum, 99.8% [95% C.I. 99.3%, 100%] for plasma and 99.7% [95% C.I. 99.0%, 100%] for whole blood (Trinity Biotech 2003).

In 968 samples from individuals of unknown HIV-1 status from a high risk population (level 3b evidence) who were subsequently found to be HIV-1 sero-negative by ELISA, two were reactive by plasma and whole blood and three by serum when tested with Uni-Gold™ Recombigen® HIV. Specificity in a high risk population was 99.7% [95% C.I. 99.0%, 99.9%] for serum, 99.8% [95% C.I. 99%, 100%] for plasma and 99.8% [95% C.I. 99 %, 100%] for whole blood (Trinity Biotech 2003).

In 2001 there were a total of 778 new cases of HIV infection in Australia (Communicable Diseases Network Australia - National Notifiable Diseases Surveillance System, 2004). Although the incidence rate in Australia is low the use of rapid of the test may be advantageous in hospital settings for quick assessment of emergency cases and in worker needle stick injuries.

CONCLUSION:

There is limited, level 2b evidence available from a single source (manufacturer). In addition, there is likely to be limited uptake of the Uni-Gold™ Recombigen® HIV test in Australia.

HEALTHPACT ACTION:

It is therefore recommended that this technology be monitored.

SOURCES OF FURTHER INFORMATION:

Constantine, N. T., Sill, A. M. et al (2003). 'Improved classification of recent HIV-1 infection by employing a two-stage sensitive/less-sensitive test strategy', *J Acquir Immune Defic Syndr*, 32 (1), 94-103.

Trinity Biotech 'Summary of Safety and Effectiveness'. [Internet]. Available from: <http://www.fda.gov/cber/pmasumm/P0300250S.pdf> [Accessed 25th February 2004]

Phillips, B., Ball, C. et al (2001). *Oxford Centre for Evidence-Based Medicine Levels of Evidence (May 2001)* [Internet]. Centre for Evidence-Based Medicine, Oxford, UK. Available from: http://www.cebm.net/levels_of_evidence.asp [Accessed 28th January 2004].

SEARCH CRITERIA TO BE USED:

HIV Infections/ classification/ diagnosis
HIV Seropositivity/ diagnosis
Immunoenzyme Techniques
Sensitivity and Specificity